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EW3039 Wrist Blood Pressure Monitor

- This instrument is an automatic electronic sphygmomanometer for measuring wrist blood pressure. The values of blood pressure and pulse can be easily measured by using oscillography.
- Applicable to the measurement of human body's blood pressure and pulse without traumatism.
- The wrist circumference that can be measured with this Wrist Blood Pressure Monitor is about 12.5cm-22.0 cm.

Our promise

Any malfunction occurring within of purchase will be dealt with according to the rules in Repair Card attached to the product.

Packing list

- ① Main unit (with pressure cuff) 1
 - ② Storage case 1
 - ③ AAA-size batteries 2
 - ④ User manual 1
 - ⑤ Repair card 1
 - ⑥ Certificate 1
 - ⑦ Record card for blood pressure measurement 1
- Included batteries are for device testing. Please verify the functions of the device during purchase.

Content

Important Safety Instructions	31
Parts Identification	32
Preparation	33
• Preparations Before Measurement	33
• Setting the date and time	34
• Positioning The Cuff	35
How to Use	36
• Instructions for correct measurement	36~37
• Methods of measurement	38~41
• Save and read measured values	42~44
• Comparison of average values of highest morning and evening blood pressures.	45
• Set (confirm) prompt time for measurement.	46~48
Other	49
• After Use	49
• Cautions During Measurements	50
• Q&A	51~53
• Care And Maintenance	54
• Troubleshooting	55~56
• Specifications	57

IMPORTANT SAFETY INSTRUCTIONS

- * Before using this instrument, please read the "IMPORTANT SAFETY INSTRUCTIONS" and be sure to use the instrument properly.
- * The cautions in this section are intended to instruct you in the proper use of this instrument and prevent danger or injury to yourself or other people.

All of the cautions contain important information regarding safety and must be followed.

WARNING: Failure to follow the caution may result in death or serious injury of people.

CAUTION: Failure to follow the caution may result in injury to people, or in damage to property.

WARNING:

Persons with serious obstructions of blood flow in the arm should consult a doctor before using this instrument. Failure to do so may result in deterioration of physical condition.


CAUTION:

- Do not make judgements yourself regarding the measurement results or self-treat yourself.
- There will be large differences between wrist blood pressure values and upper-arm blood pressure values for persons with periphery circulatory obstruction due to diabetes, liver disease, kidney disease, arteriosclerosis, high blood pressure, etc.
- Take measurements according to your doctor's instructions and have your doctor do the diagnosis.
- Take medicine according to your doctor's instructions.
- Do not use for purposes other than the measurement of blood pressure. Doing so may result in an accident or malfunction.
- Do not use this instrument on infants or on people who cannot express their feelings. Doing so may result in problems or an accident.
- Do not disassemble this instrument, or attempt to repair or modify it yourself. Doing so may result in an accident or malfunction.

SAVE THESE INSTRUCTIONS

PARTS IDENTIFICATION

Please check the instrument.


<About batteries>

- Panasonic alkaline batteries (two AAA-size batteries) can provide approximately 300 uses (for measurements 3 times a day, approximately 3 months of use). (Measurement conditions: Room temperature: 23°C Pressurized to 170mmHg; Wrist circumference: 17cm)
- Batteries are included only for testing purposes and with no guarantee of battery life. (In some areas batteries are not included with the unit.)
 - If the measurement notification sound is used, battery life will be shorter.
 - At lower room temperatures, battery life will be reduced.
 - When using batteries other than alkaline batteries, the number of uses may be greatly reduced. Always use alkaline batteries.
 - When installing batteries, both batteries should be of the same type and from the same manufacturer, and should be installed together at the same time.
 - If the instrument will not be used for a long period of time, remove the batteries. Leaving the batteries installed may result in battery leakage and damage to the instrument.
 - Do not discard batteries with normal garbage. Take to recovery center.

PREPARATIONS BEFORE MEASUREMENT

* Always use alkaline batteries.


1 Remove the unit from its storage case.

Inserting dry cell batteries.

- (1) Remove the battery case by pushing it in the direction of the arrow.
- (2) Insert the batteries, ensuring the correct polarity.

2 Slide the battery case back into place on the main body until it clicks into place.

Replace batteries when:


- The mark shown at right appears during measurement.
- No display is shown when the Start/Stop button is pressed.

⚠ CAUTION:

- Always follow the cautions printed on the batteries.
 - Check the polarities of the batteries and install them correctly.
 - Promptly remove exhausted batteries.
 - If the product will not be used for a long time, remove the batteries.
 - Batteries should be used before their expiration date.
- Use after the expiration date may result in injury or staining of the surrounding area due to generation of heat, bursting, or leakage.

SETTING THE DATE AND TIME

* Please be sure to set the date and time before taking the first measurement after purchasing the blood pressure monitor or after changing the batteries. (If a measurement is taken without setting the date and time, r:--:-- will be displayed.)

1 Hold the **設定** button on bottom surface of the monitor pressed **for at least 3 seconds**. The day value will blink.

If the date and time have already been set and the button is not held pressed long enough, the display will change to the measurement notification sound setting display.

2 Set the day.

Press the **変更** button repeatedly to set the correct day.

3 Set the month.

Press the **設定** button. The day value will blink. Press the **変更** button repeatedly to set the correct month.

4 Set the hour.

Press the **設定** button. The hour value will blink. Press the **変更** button to set the correct hour.

5 Set the minute.

Press the **設定** button. The minute value will blink. Press the **変更** button repeatedly to set the correct minute.





6 Press the **設定** button at the exact time the minute changes **to finalize** the date/time setting.

7 Press the **开始/结束** button to **switch off** the power.

When date and time are blinking, the clock has stopped.

Feb 29 of leap year is displayed as Mar 1, so please reset the date. Please use the tip of a ballpoint pen to press the buttons of "設定" and "変更" lightly.

POSITIONING THE CUFF

1 Attach on the inside of the wrist.

- Wrap the cuff on bare skin.
- Measurable wrist circumference: 12.5cm to 22.0cm

2 Fasten the cuff firmly

- If the cuff is wrapped loosely, the measured value for blood pressure will be high.
- Avoid twisting the belt when fastening.

When measuring on the right wrist:

Position the blood pressure monitor as shown in the figure.

The blood pressure monitor can be used by wrapping the cuff on either wrist. However, since it has been said that the difference in blood pressure between the left and right wrists may be around 10mmHg, be sure to always use the same wrist for measurements.


FOR ACCURATE MEASUREMENTS

Take measurements at around the same time every day, using the same wrist and in the same posture.

- When taking measurements in the morning, take the measurement immediately after waking up when there is no urge to urinate, or if that is difficult, take measurements before eating breakfast without moving your body.


Before measurements, be sure to sit quietly for 4 to 5 minutes.

When measuring while sitting down:

1 Rest your elbows on a table.

2 Position yourself so that the blood pressure monitor is at the same level as your heart.


3 With the palm of the hand facing up, release any tenseness in your muscles.

If measurements are taken with the wrist or hand in the following positions, accurate measurements will not be obtained.

Wrist is bent inwards.

Wrist is tense and fist is being made.

FOR ACCURATE MEASUREMENTS(CONTINUED)

If the blood pressure monitor is positioned lower than the heart, the measured blood pressure value may be high.

If the posture applies pressure to the stomach (leaning forward), accurate measurements may not be possible.

Do not move your body, your arm, or your fingers while taking measurements.

Do not use the blood pressure monitor near cellular phones or other devices which emit electromagnetic waves. Doing so may result in malfunction.

If the blood pressure monitor is positioned higher than the heart, the measured blood pressure value may be low.

Do not touch the blood pressure monitor while taking measurements.

Do not talk while taking measurements.

TAKING A READING

Stay quiet and relax while taking measurements.

1 Place the storage case under your arm. (Refer to "FOR ACCURATE MEASUREMENTS" on see pp. 36 and 37.)

2 Press the 开始 / 结束 button.

- All displays will be shown for approximately 2seconds.
- Pressurization will start automatically.
- The heart mark ♥ will blink with sound when pulse is detected during pressurization.
- If the measurement notification time is set, the bell mark will be shown. (To set or check the measurement notification time, see pp. 46 to 48.)

Since measurements are performed during pressurization, there may be changes in operating sounds or the pressurization rate during pressurization. This is normal and is not a malfunction.

3 When measurement has been completed, the beeps will sound and **the blood pressure and pulse readings will be displayed.**

- Pressure will be automatically released from the pressure cuff.
- The blood pressure indicator lamps will light (for approximately 6 seconds).
- If the measured value is in the high blood pressure region, the values will blink (for approximately 6 seconds). (For how to read the blood pressure indicator lamps, see p. 41.)
- The \bar{H} or \bar{L} mark will be lit up. (For measurement values which can be stored, see p.42~44)

* The pulse rate value is the beats per minute calculated based on the pulse rate during measurement.

4 To switch off the blood pressure monitor, **press** the 开始/结束 button.

- If you forget to switch off the blood pressure monitor it will switch off automatically after approximately 5minutes.

TAKING A READING (CONTINUED)

MOVEMENT DETECTION FUNCTION

If the wrist bound with sleeve strap flexes unconsciously during the measurement, the blood pressure results may be inaccurate at times.

- The "Movement Detection Function" is a function that detects excess pressure on the wrist strap caused by movement of the user's hand or wrist during blood pressure measurement and shows the "Movement Detection" symbol in the display when this occurs.
- If the "Movement Detection" symbol blinks, please resume the proper position of your arm (see pp.36 and 37).
- The "Movement Detection" symbol will blink if the body moves during blood pressure measurement. When the result is displayed and body was moving during measurement, the "Movement Detection" symbol will also blink.

(When a saved result, measured while the body was moving, is selected, the "Movement Detection" symbol will also blink.)


"Movement Detection" symbol display	Alarm	Explanation
	Long "peep-peep" sound	Body movement detected <ul style="list-style-type: none"> • Involuntary flexure of hand or wrist. • Overexertion of the hand or wrist. • Movement of the hand or wrist etc. Please resume the right position and measure again(see pp.36 and 37).
	Short "peep-peep-peep" sound	Movement that may cause inaccurate measurement result occurred: <ul style="list-style-type: none"> • Large degree of movement • Multiple movements • Large degree of movement of the hand or wrist etc. When "E" is displayed on the screen, please measure again. ("E" will be displayed when:) <ul style="list-style-type: none"> • The fastener of sleeve strap is loose.

*This is a reference function that designed for accurate measurement. In order to obtain more accurate blood pressure measurement results, please measure for 2-3 times under static status even if the "Movement Detection" symbol does not blink (see p.36).

IRREGULAR PULSE DETECTION FUNCTION

If the heart rate measured exceeds the average value for more than ± 25% during measurement, the symbol of "Irregular Pulse" will be displayed to remind that attention is needed after measurement.

- After measurement, the symbol of "Irregular Pulse" will be displayed after 6 seconds of blinking.
- When the symbol of "Irregular Pulse" is displayed, please measure again.
- If the heart rate changes too much during measurement, it may affect the accuracy of measured value.

* Blood pressure values are displayed alternately in mmHg and kPa units.

- When you save and select the value measured with irregular pulse, the symbol of "Irregular Pulse" will be displayed as well.

This symbol is not designed to indicate irregular pulse.



- This is a reference function intended for ensuring the accuracy of measured value. In order to ensure the accuracy of measurement, please measure for 2-3 times under static status even if the symbol of "Irregular Pulse" is not displayed.(see.p.27)
- If the symbol appears constantly, you should consult a doctor regarding your health condition.

Please do not estimate the result and receive treatment by yourself: doctor's advice must be followed.

TAKING A READING (CONTINUED)

Reading the blood pressure display

Whether the systolic and diastolic blood pressure values are in the high blood pressure region or normal blood pressure region are judged based on the WHO/ISH(Word Health Organization/International Society of Hypertension)blood pressure classifications,and the blood pressure value blinks if required.

If the measured blood pressure value is in the high blood pressure region, the values will blink for approximately 6 seconds.

- Systolic blood pressure:18.7kPa (140mmHg) or higher
- Diastolic blood pressure:12.0kPa (90mmHg) or higher (Only immediately after measurement) (1kPa=7.5mmHg)

If measurements could not be taken properly

The remeasure mark "E" will be displayed.

When taking successive measurements,press the 开始/结束 button to switch off the blood pressure monitor once and stay quiet for 4 to 5 minutes,and then take the next measurement.

Advice:

Often with the second measurement,there is less tension and you feel more relaxed,so the blood pressure value is lower. For people who easily get tense,measuring twice each time is very effective. When two measurement are taken,normally the smaller measurement is stored.

When taking two measurements,always stay quiet for 4 to 5 minutes after the first measurement before taking the second measurement.



※1 1999 World Health Organization-International Society of Hypertension Guidelines for the Management of Hypertension

STORING AND RECALLING MEASUREMENT VALUES

Storing measurement values

Up to 90 readings each can be stored in the memory for two people. Stored data will remain in memory even if batteries become exhausted or batteries are removed.

If measurements were taken without setting the date and time, the date and time will be shown as "-/- -:-".

1 After measurement has been completed, the mark will be lit. Press the button for the memory in which you want to save the reading.

E.g. A....Recorded in
B....Recorded in

2 Storing is completed.

- Up to 90 measurement values can be stored. When the memory is full and additional values are stored, the oldest value in memory will be deleted and the new value will be stored. At that time, the data numbers for the measurement values will change.

- If measurements could not be taken properly (the mark is shown instead of the blood pressure value), measurement values cannot be stored. Data numbers will not be changed.

* However, if only the pulse rate value shows the , measured blood pressure values can be stored.

- The memory number, measurement date, and time are alternately displayed approximately every two seconds.

STORING AND RECALLING MEASUREMENT VALUES (CONTINUED)

Recalling measurements


Measurement values stored in memory can be recalled at a later time or after the power has been switched off.

1 Press the 记忆/调出 button of the person for whom you wish to recall a reading.

- The average of all measured values in memory will be shown.
- The mark will be lit.
- For example, if 90 measured values are stored, the average of 90 measurements will be shown. (If only 1 measurement is stored, the mark will not be shown.)
- Blood pressure indicator lamp will be lit up for approximately 6 seconds.

2 If you press the 记忆/调出 button again, the most recent stored measurement value will be shown.

- Each time the 记忆/调出 button is pressed, the stored measurement values will be recalled in order, starting with the most recent.
- Memory number, date and time of the reading will be displayed for approximately two seconds each.

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3 To switch off the blood pressure monitor, **press** the 开始 ▶ 结束 button.

If you forget to switch off the blood pressure monitor, it will switch off automatically after approximately 30 seconds.

To erase all measured values stored in memory:

Press the 记忆/调出 button to recall stored measurement values.

Press the 记忆/调出 button again and hold it pressed until the measurement values disappear (longer than approximately 3 seconds).

* It is not possible to erase single measurement values.


When all measurement values have been deleted, the selected mark or mark and mark will be displayed (with sound).

COMPARISON OF THE AVERAGE VALUES OF THE HIGHEST MORNING AND EVENING BLOOD PRESSURES.

This function displays the average values of the highest morning and evening blood pressures, and the measured values of different periods of morning and evening can be compared.

This function Comparing morning and evening values can be helpful to detect morning fast-increasing pressure type of hypertension as soon as possible.

However, the comparison will be unavailable if 日期 / 时间 is not set.

1 Press the 早晚比较 button.

- Display the average values of the highest morning blood pressure and the highest evening blood pressure of User 1.

Morning : 3:00 AM-11:00 AM
 Evening : 18:00 PM - 2:59 AM

- Display saved average value of the highest blood pressure during above periods.
- If time is not set or saved values are less than twice, the measured value is displayed as [- -].

2 Press the 早晚比较 button again.


- Display the average values of the highest morning and evening blood pressure of User 2
- Every time the button is pressed, the average values of the highest morning and evening blood pressures for User 1 and User 2 are displayed alternatively.

3 When measurement is finished, press 开始▶结束 buttons.

- Even if you forget to shut off power supply, it will be shut off automatically after 30 seconds.

SETTING (CHECKING) THE MEASUREMENT NOTIFICATION TIME

To let you take measurements at the same time every day, 3 times per day can be set as notification times.

1 Press the 設定 button on the bottom surface of the blood pressure monitor.

- Notification Time 1 will blink.

* If date and time have not been set, or the button is held pressed for more than 3 seconds, the display will change to the 日期/时间 setting display.

2 Press the Set button again.

- The \triangle and OFF will blink.

3 Press the Change button to display ON.

- The \triangle and ON will blink.
 - ON \rightarrow Notification sound will sound.
 - OFF \rightarrow Notification sound will not sound.
- When ON is displayed and the Set button is pressed, the time becomes 0:00 and the \triangle and ON will light.

- ### 4
- ① Press the Set button. The hour value will blink.
 - ② Press the Change button repeatedly to set the desired hour.

SETTING (CHECKING) THE MEASUREMENT NOTIFICATION TIME (CONTINUED)

- 5**
- ① Press the Set button. The minute value will blink.
 - ② Press the Change button repeatedly to set the desired minute.
- To set only Notification Time 1, press the 开始 / 结束 button to complete the setting.
 - To set Notification Times 2 and 3, follow steps 6 and 7 below.

- 6** Press the Set button. Notification Time 2 will blink. Repeat steps 2 to 5 to set the desired notification time.

- 7** When setting of Notification Time 2 has been completed, press the Set button. Notification Time 3 will blink. Repeat steps 2 to 5 to set the desired notification time.

- 8** When setting of Notification Time 3 has been completed, press the Set button. The times set for
- 9** To complete the setting, press the 开始 / 结束 button. Notification Times 1 to 3 will be shown automatically at 2 second inter vals.

To change the notification time setting:

- 1** Press the Set button on the bottom surface of the blood pressure monitor. (Notification Time 1 will blink.)
 - To change the setting for Notification Time 2 or 3, press the Change button to display Notification Time 2 or 3.
- 2** Press the Set button twice. The hour value will blink.
- 3** Press the Change button repeatedly to set the desired hour.
- 4** Press the Set button. The minute value will blink. Press the Change button repeatedly to set the desired minute.
- 5** Press the 开始/结束 button to complete the setting.

To change the measurement notification sound:

- 1** Press the Set button on the bottom surface of the blood pressure monitor. (Notification Time 1 will blink.)
 - To cancel notification for Notification Time 2 or 3, press the Change button to display Notification Time 2 or 3.
 - 2** Press the Set button once. The  and  will blink.
 - 3** Press the Change button.  will be displayed.
 - 4** Press the 开始/结束 button to complete the setting.
- When the set time is reached, the measurement notification sound will sound to notify you that it's time to take a measurement.
 - If the set time is reached in the middle of taking a measurement, the measurement notification sound will not sound.
 - When the measurement notification sound sounds, it can be stopped by pressing the 开始/结束 button. If the button is not pressed, the sound stops after 1 minute.
 - If the blood pressure monitor is put into a bag, desk, etc., it may not be possible to hear the measurement notification sound.


Advice:

If a record table is made of the measurement values taken around the same time, you can easily see your own daily blood pressure changes.

AFTER USE


Storage

1

1 Fold the pressure cuff as indicated.

2


2 Put the blood pressure monitor in the storage case.

Avoid pinching the cuff with the storage case.

CAUTIONS DURING MEASUREMENTS

In the following situations, accurate measurements cannot be taken.

Within 1 hour of eating or drinking alcohol

Immediately after exercising or taking a bath

Wait at least 20 minutes and take a measurement while staying quiet.

Immediately after drinking coffee or tea, or smoking

In cold places

Measurements should be taken at room temperatures of around 20°C.

In a moving vehicle

When you need to urinate

Wait several minutes after urinating before taking measurements.

Difference between wrist and upper-arm blood pressure values

There may be large differences between wrist blood pressure values and upper-arm blood pressure values for persons with periphery circulatory obstruction due to diabetes, liver disease, kidney disease, arteriosclerosis, high blood pressure, etc.*1

Because of this, do not make judgements regarding blood pressure yourself. Always follow your doctor's directions.

*1 Persons with basic diseases and poor circulation in the wrist.

Q & A

Q Why are the blood pressure values measured at home different from those measured at the doctor's office?


A Blood pressure varies 24 hours a day. In addition, it can change greatly depending on the weather, your mental state, exercise, etc.

In particular, there is the phenomenon of "White uniform hypertension" at doctor's office, where the blood pressure measured in a doctor's office is higher because of stress or worry.

In addition, if the blood pressure measured at home is higher than that measured at the doctor's office, the following points should be considered and care should be observed where necessary:

- ① The blood pressure monitor is not at the same height as the heart.
If measurements are taken with the blood pressure monitor lower than the heart, the measured blood pressure values will be higher.

- ② The pressure cuff is wrapped loosely.
Wrap the pressure cuff securely around the wrist.


- ③ Do you feel worried or irritated?
Before taking measurements, take 2 or 3 deep breaths and relax to stabilize blood pressure. In addition, it is recommended that you sit quietly for about 4 to 5 minutes before taking measurements.

Q & A

Q Why are measured blood pressure values different each time I take a measurement?


A

- ① Blood pressure varies throughout the day, and can also vary depending on how the pressure cuff was wrapped or your posture. Because of this, measurements should be taken at the same time and under the same conditions every day.
- ② For persons taking blood pressure medication, blood pressure may change greatly due to the effects of the medication.
- ③ When taking 2 measurements in succession, be sure to stay quiet for 4 to 5 minutes after the first measurement before taking the next measurement.

Q Why are the measured values different for the left and right wrist?

A

It is said that differences in blood pressure values between the left and right wrists may be around 10 mmHg. The blood pressure monitor can take measurements on either the left or right wrist, but all measurements should be taken using the same wrist.

Q The "Movement Detection" symbol is not displayed while the hand or the wrist is moving. why?

A The "Movement Detection Function" is a function that detects excess pressure on the sleeve strap which is caused by movement of the user's wrist during measurement and displays the "Movement Detection" symbol when this occurs. Because of this, the "Movement Detection" symbol may not be displayed if the movement does not affect the sleeve strap.

If you are unsure about the measurement, please do it again.

Q & A

Q The measured value is unusually higher although the position of the hand and the wrist is fixed according to the "Movement Detection" symbol. Why?

A The "Movement Detection" symbol is just a reference indication for correct measurement posture (see pp.39).

If the measured value is unusually higher or lower after the symbol of "Movement Detection" is displayed, please measure again.

If the above situations occurs, the following factors may exist:

- Incorrect measurement method (the binding method of sleeve strap, posture, etc)?

Please refer to the methods of pp.36 ~ 37.

- Measurement immediately after smoking.
- Measurement when user's mood is not steady.
- Measurement in cold circumstance.
- The change of user's pulse is small.

If you are unsure about the measured value, please measure again after 4-5 minutes rest.

Q When should I measure my blood pressure?

A We recommend that you measure your blood pressure 3 times a day.

First time: When you get up in the morning. (This is when your blood pressure is most typical of its value during the night and when the movements of your body have had hardly any effect on it.)

Second time: Before eating lunch. (This is when your blood pressure is not affected by food and when it is most typical of the highest level of the day.)

Third time: When you go to bed at night. (This data value is roughly typical of the average blood pressure value for the day.)

Get to know your blood pressure pattern by measuring your blood pressure in the morning, at midday and in the late evening, and use the data you obtain to manage your blood pressure.

Suggestions of use methods (Please refer to page 45~ 47)

- Please set the prompt time as the time mentioned above when using.
- By using the function of comparing morning and evening values, it is helpful to detect morning fast-increasing pressure type of hypertension as soon as possible.

CARE AND MAINTENANCE

Do not disassemble the instrument or attempt to repair or modify it.

Always wrap the pressurization cuff around the wrist for use. Do not bend the pressurization cuff backwards.

- Doing so may damage the pressurization cuff.

Do not apply excessive force to the instrument or drop the instrument.

- Doing so may cause damage.

Do not insert dust or foreign objects into the instrument.

- Doing so may result in damage.

If the instrument is stored at temperatures below the freezing point, do not use it immediately. Leave it in a warm place for at least an hour before use.

- If the instrument is not allowed to warm up, it may not pressurize.

Do not touch the instrument during use.

- Doing so may prevent accurate measurements.

If the instrument becomes dirty, clean it with a soft cloth moistened with warm water or soapy water. (Do not use alcohol, benzine, or paint thinner.)

- Use of such chemicals may result in cracking or discoloration.

If the instrument will not be used for a long period of time (more than 1 month), remove the batteries.


- Leaving the batteries installed may result in battery leakage and damage to the instrument.

The pressure cuff cannot be washed.

For storage, avoid high temperatures, high humidity, and direct sunlight.

- Storage under such conditions may result in damage.

TROUBLESHOOTING

Display	Status before error	Cause and solution
Remeasure  mark is displayed. 	Pressurization to more than 280 mmHg was performed.	Were measurements taken according to the proper procedure and in the correct posture? See pp. 36 ~ 37.
	Pressure decreased suddenly and the remeasure mark  appeared.	Was the pressurization cuff wrapped properly? See p. 35.
	The heart mark  blinked only a few times or not at all.	
	No pressurization	
	Your hand or wrist moved.	Has the unit detected movement? See p. 39.

Symptom	Possible cause
Systolic or diastolic blood pressure value is high or low.	<ul style="list-style-type: none"> • During measurement, the level of the instrument was higher or lower than the level of the heart. • The pressure cuff was not wrapped properly. • You moved or talked during measurement. See pp.35 to 37.
Value is different from that measured by the doctor.	<ul style="list-style-type: none"> • Blood pressure varies in response to minute changes in your mental state, such as your reaction to the doctor's presence.
Measurement values are different with each measurement.	<ul style="list-style-type: none"> • Relax 4 or 5 minutes, and take another measurement. See pp. 36 ~ 37.

Symptom	Possible cause
No display appears even when the 开始 / 结束 button is pressed.	<ul style="list-style-type: none"> • Batteries are installed backwards. See p. 33.
When the 开始 / 结束 button is pressed, the display disappears during pressurization ( mark appears).	
During pressurization, the sound or pressurization rate changes.	<ul style="list-style-type: none"> • Since measurements are performed during pressurization, there may be changes in operating sounds or the pressurization rate during pressurization. This is normal and is not a malfunction.
The  mark flashes when readings are being taken.	<ul style="list-style-type: none"> • There was a significant fluctuation in pulse when readings were being taken. See p. 40.
The  mark flashes when readings are being taken.	<ul style="list-style-type: none"> • Your hand or wrist moved. See p.39.
The  mark flashes even though I did not move my hand or wrist.	<ul style="list-style-type: none"> • In some cases this will be displayed when you have tensed the muscles in your arm. Relax the muscles in your arm and take the reading again. See pp. 36 ~ 37.
Damage of the main part	<ul style="list-style-type: none"> • Do not repair it yourself. Please contact our designated service agency for wrist blood pressure monitors to examine and repair the device.

***The warranty of this product is only valid in the PRC (Not available in Hongkong, Macao and Taiwan for the time being)**

Specification

Power	DC3V AAA-size alkaline batteries × 2
Display mode	Digital LCD
Measurement mode	Oscillometric system
Measurement range	Pressure: 0kPa~37.3kPa(0mmHg~280mmHg) Pulse rate: (30~160)beats/minute
Accuracy	Pressure: 15℃~25℃ within ±0.4kPa(3mmHg) 5℃~40℃(out of 15℃~25℃) within ±0.8kPa(6mmHg) Pulse rate: within ±5%
Rate of pressurization	0.5kPa/s~1.1kPa/s(3.75mmHg/s~8.25mmHg/s)
Operating condition	Temperature 5℃~40℃ Relative humidity ≤ 80%
Storage & transportation condition	Temperature -20℃~55℃ Relative humidity 10% ~ 93%
Measurable wrist circumference	About 12.5cm~22.0cm
Weight	Approx. 116g (Excluding the batteries)
Dimensions (Main part)	Approx. 78mm × 63mm × 26mm
Protection against electric shock	Internally powered equipment, Type B applied part
Pressurization mode	Pressure-regulated mode of oscillating-mode pump
Fast air-discharging mode	Fast-open air-discharging valve mode
Pressure detection	Semiconductor pressure sensor
Pulse detection	Semiconductor pressure sensor

Certificate: